



Top 5 Benefits of MobiLab® Desktop

MobiLab Desktop is a barcode specimen collection solution that runs on both PC workstations and laptops on mobile carts. Our customers recently shared with us the top five reasons they use MobiLab Desktop in addition to their MobiLab handheld devices:

- 1. Eliminate Mislabeled Specimens.** By using hardware that is already part of their workflow and intuitively designed software to print labels at the patient bedside, teams can eliminate mislabeled specimens associated with the manual workflow errors of the past. No more label batch, no more walking to get labels from a distant printer, no more mislabeling of specimens.
- 2. Less Hardware Costs.** Phlebotomists find using MobiLab on handheld devices perfectly fits their needs. MobiLab Desktop is used from a computer desktop, which means no additional hardware is needed for nurses and other caregivers who primarily work on a PC, helping your team save on hardware costs across the board.
- 3. Training & Adoption.** An intuitive design means decreased training time. Having the application available on the same device already used for work improves adoption for specimen collection. Caregivers who primarily use a PC in a patient room or a workstation on wheels (WOW) for documentation, prefer a solution that runs on the same device they already use for their workflow.
- 4. Outpatient Areas.** MobiLab Desktop is a good fit for outpatient areas that already have a PC workstation next to each draw chair. This allows nurses and others doing collections to remain at the patient's side, rather than moving constantly to complete tasks in their workflow.
- 5. Specimen Collection Tracker.** MobiLab Desktop supports a view-only display of the draw list, allowing staff to track pending collections and those that have not yet arrived in the laboratory. This also lets leaders easily see if additional support is needed for collection and allows techs to prepare for incoming work.

Start achieving these benefits today. Schedule a demonstration of MobiLab Desktop by contacting us using the information below.